

Humid air leakage into enclosure after N₂ shutoff at port #3

