


BPM Performance at AWAKE






Spencer Gessner
12 December, 2018

21 p⁺ BPMs in TT41



16 "upstream" BPMs monitored by SPS OP

5 "downstream" BPMs monitored by AWAKE


From TT40/TT41/SPS




Plasma Cell



BPM 412425
After Plasma Cell



Detail B - Scale 1:25



Detail A - Scale 1:25

NOTA :
Layout from Plasma Cell to Shielding Wall (final part of Beam Line), see drawing SP5LAWK_0009

Layout of upstream part of TT41 Beam Line see drawings SP5JTT4100003 and SP5JTT4100004

Layout of Electron Beam Line, see drawing SP5LAWK_0032

DES/DR: F. GALLEAZZI
 CONTROLLED BY: Y. MUTTI
 RELEASED BY: A. PARKING
 APPROVED BY: F. GALLEAZZI
 CAD Document Number: ST0644812_02
 REPLACES: SP5LAWK_0008
 DATE: 2015-03-27
 NOM/NAME: F. GALLEAZZI
 ZONE: MAJ Ligne Proton + Ajout Elements Vide
 MODIFICATION: Correction Longueur BTV - BCTF - BPM

IND.	DATE	NOM/NAME	ZONE	MODIFICATION
D	2017-06-14	F. GALLEAZZI		MAJ Ligne Faisceau - Version Juin 2017
C	2015-07-10	F. GALLEAZZI		MAJ suite commentaires lors de l'approbation
B	2015-06-01	F. GALLEAZZI		Correction Longueur BTV - BCTF - BPM
A	2015-03-27	F. GALLEAZZI		MAJ Ligne Proton + Ajout Elements Vide

Machine and installation layouts for AWAKE		DES/DR: F. GALLEAZZI	2014-11-28
TT41 - AWAKE		CONTROLLED BY: Y. MUTTI	2017-07-27
LAYOUT QTF 412000 - PLASMA CELL		RELEASED BY: A. PARKING	2017-07-28
TT41 - AWAKE		APPROVED BY: F. GALLEAZZI	
LAYOUT QTF 412000 - PLASMA CELL		CAD Document Number: ST0644812_02	
PROJECT ENGINEER		FOR INFORMATION	REPLACES: SP5LAWK_0008
		SCALE: 1:50	IND. 1


Proton BPMs

No issues with software. No issues with hardware or readout for “upstream” BPMs.

“Downstream” BPMs give “bad” readings in the following conditions:

1. High density Rb vapor in the line.
2. Electron beam present.
3. Proton pre-bunches.

“Bad” BPM readings




Bad readings are obvious, unphysical deviations of the beam trajectory. All 5 downstream BPMs are affected, but the final 2 are most strongly affected.

Proton Pre-bunches

In July, Livio observed significant charge ahead of the main pulse.


This also led to bad BPM readings. We alerted SPS OP and they changed RF voltage and kicker timing to avoid this problem.

Manfred also adjusted the self-trigger threshold for the BPMs.






Was this an exceptional case/accident? Do we need a pre-bunch monitor?

Plasma Effect on Proton BPMs







M. Wendt

Since April 2017, we have observed the effect of the Rb vapor/plasma on the AWAKE BPMs. Manfred proposes to use stripline BPMs instead.



12 e⁻ BPMs in TT43/TT41

- 7 BPMs in TT43
- 5 BPMs in TT41 common to p⁺
- p⁺ BPMs



NOTA : DETAILED LAYOUT
1st Part (Photoinjector + Period 430000-430100)
see drawing SPSLAWK_0014
2nd Part (Period 430100-430200 - End Line)
see drawing SPSLAWK_0012

Revision and Installation Sheets for AWAKE		SCALE	DATE	BY	APPR
TT43 - AWAKE		1:10	2017-04-12	AWAKE	AWAKE
ELECTRON BEAM LINE LAYOUT					
TT43 - AWAKE					
LAYOUT LIGNE FAISCEAU ELECTRON					
PROJECT LEADER (BY INFO)		INFORMATION		SPSLAWK_0032	

IND.	DATE	NOM/NAME	ZONE	Modifs sans remarques de l'Approbation / Modif. de l'Approb.

Electron BPMs

We hoped for simultaneous electron-proton readings on these BPMs.

- Maybe this was never possible. . .


On the plus side, the readings appear to be very precise and not strongly affected by the presence of plasma.

We had issues with timing and synchronization of readout of BPM data.


- Athanasios solved this problem but I think we need to revisit triggers, time-stamping, and synchronization.

Electron-Proton BPM Correlations

X-plain, eBPM v.s. pBPM, tilt 0.685



Y-plain, eBPM v.s. pBPM, tilt 0.685



Misha, Shengli, Lars, Michal, did cross correlations on electron and proton BPMs. This was extremely useful and should be part of system validation in the future.