

Report: 24

Activity: SPS

*Chiara Pasquino,
Meeting 02/03/2020*




LS2



Activities of the past 2 weeks

- BA1:
 - Beamline installation of 135, 136 ,142.
 - Pumpdown of 111-112; → postponed;
- BA2:
 - Magnet campaign;
- BA3:
 - Installation of the 5th cavity sector; 1st cavity sectors installation completion with sector valves; Couplers installation ongoing on 3rd cavity sector;
- BA4:
 - FWS installation;
- BA6:
 - Pumpdown of 661 and 662; → postponed (linked with BA1)
- TI2 – TI8:
 - Pumpdown of TI8 line; → vented due to alignment issues on a pumping module and a BTV;

Few Pics from last meeting..




TT66 – HiRadMat
Survey smoothing
campaign..
Sector vented, flanges
disconnected, waiting for
different type of supports
from MSC to replace the
old ones. 4 magnets
involved in this
exchange.

ION PUMPS AND SAFETY..




Abel reinstalled them, BE-RF is warned about it.

Position of the ion pumps local boxes discussed in situ with integration and BE-RF... still someone removed them while being powered!



LSS1 installation




- S135 – mechanically closed, pumping down (leaks!)
- S136 – sector valve position to be modified then ready for final connection
- S142 – vacuum drifts and instrumentation are installed, waiting for final alignment of the mask to complete the mechanical installation

Total doses are OK, personal under control, dose sharing in place and we will rotate during this week too.

LIU-Projects: LSS1, LSS3 & LSS5 status

- LSS1: possibly starting pumping 135 & 136.
- LSS5: production completed, sorting in SMA18 this week + transport of last 159 to be coated.
- LSS3:
 - 6th cavity sector being installed;
 - Integration clash in two positions between sector valve piston and magnet pumping port. Solutions: changing valve type (different version of VVSBs) or swapping position of the valve. → solved with a different version of VVSB and different configuration of the ion pump.
 - Second time we find a leak on brand new ion pump feedthrough.


TDC2 tunnel inspection



TDC2 – LS2 unforeseen actions

- First Collimator exchange;
- Exchange of the vacuum chambers downstream both sets of splitters due to corrosion (≈ 10 chambers).
- Needs an ECR cause some chambers have been welded in place at their first installation.
- Sector with Splitters 2 will have to be vented (risk of opening existing varnished leaks).

DASHBOARDS - 2020



Activities of the next 2 weeks

- BA1:
 - Pumpdown of 111-112;
 - Pumpdown of 135-136-142;
- BA3:
 - Installation of the 6th cavity sector; 1st cavity sectors installation completion with sector valves;
- BA4:
 - FWS leak detection;
- BA6:
 - Pumpdown of 661 and 662;
- TI2 – TI8:
 - Pumpdown of TI8 line;

Resources distribution

- LSS1 reinstallation: Tony
- LSS3 reinstallation: Jarmo
- Arcs pumpdown and leak detection: Jarmo, Anthony, BINP & 40/30;
- Support to aC coating: Tony & Jarmo;
- Acceptance test & lab activities: Tony;
- Jobs follow-up & ECR update (Design & Production): Chiara & Tony;
- Stock check and purchasing: Chiara & Tony.

Thank you !




LS2




Pumpdowns

LS2 - arcs pumpdowns





LS2 - DS pumpdowns



Lengthy pumpdowns (ion pump flashing) – arcs/DS are behaving all similarly.





Pumpdowns



Lengthy pumpdowns in LS1 too, very similar pumpdown considering we are comparing 1 arc with an half arc (double pumping speed).
For info: 510 – LS1 1 year exposure;
411 – LS2 3 months exposure;

Longer exposure --> difficult flashing.

Pumpdowns – aC coated Vs uncoated sectors (once on ion pumps)



Pumpdowns - summary

- Pumpdown on TMPs is strongly dependent on:
 - Exposure time to air:
 - Location of the openings, length of the intervention;
 - It strongly affects the behaviour of the ion pumps during flashing; → is becoming a limiting factor (all SPS pumping groups are in use at the minute).
 - New ion pumps shows a more 'stable' behaviour to flashing, but still the process is lengthy in time.
- For long exposures (LS1 and LS2), pumpdowns are very similar for coated and uncoated sectors;
- Pressure recovery once on ion pumps, is very similar for coated and uncoated sectors, after long exposure.
- After LS2 I would test ion pumps performance at different exposures to air/N₂/noble gases.