

Frequency combs for differential spectroscopy of atomic clocks

Monday, 16 October 2023 12:00 (30 minutes)

Over the past 20 years optical frequency combs [1], with atomic clocks [2], have been a powerful and enabling technology in the context of time and frequency measurement [1,2]. Impressively, optical atomic clocks have yielded an 8 order of magnitude improvement in accuracy in the past 30 years. These improvements are fueling a push toward redefinition of the SI second to optical atomic references [3], as well as application of atomic clocks to tests of fundamental physics [4] and as relativistic gravitational sensors [5-6]. Unfortunately, the long measurement times needed to average down clock quantum projection noise and local oscillator noise to reach measurement stabilities at and beyond the 10⁻¹⁸ level, limit the feasibility of next-generation applications.

I will present the improved instability results for an inter-species optical atomic clock comparison using a differential measurement technique, Figure 1. In this technique, the single ion ²⁷Al⁺ clock near and the ¹⁷¹Yb lattice clock shared a common local oscillator using the phase coherent wavelength conversion with an optical frequency comb. This technique enabled nearly a factor of 10 improvement in 1-s measurement resolution and a 100-time improvement in averaging time to reach a measurement instability of 10⁻¹⁸. Improvements in the measurement stability was achieved via a minimization of laser noise aliasing, and via improvement in the ²⁷Al⁺ clock quantum projection noise by increasing its probe time by mitigating laser-atomic decoherence [7].

References

- [1] Fortier, T.M. and Baumann, E., “20 years of developments in optical frequency comb technology and applications,” *Communications Physics* 2, Article number 153 (2019).
- [2] Ludlow, A. D., Boyd, M. M., Ye, J., Peik, E. & Schmidt, P. O., “Optical atomic clocks,” *Rev. Mod. Phys.* 87, 637–701 (2015).
- [3] Riehle, F., Gill, P., Arias, F. & Robertsson, L. “The CIPM list of recommended frequency standard values: guidelines and procedures. *Metrologia*,” 55, 188–200 (2018).
- [4] BACON collaboration: “Frequency ratio measurements at 18-digit accuracy using an optical clock network,” *Nature* 591 (2021).
- [5] Rosenband, T. et al. “Frequency ratio of Al⁺ and Hg⁺ single-ion optical clocks; metrology at the 17th decimal place,” *Science* 319, 1808–1812 (2008).
- [6] Mehlstäubler, T. E., Grosche, G., Lisdat, C., Schmidt, P. O. & Denker, H. Atomic clocks for geodesy. *Rep. Prog. Phys.* 81, 064401 (2018).
- [7] Kim, M et al., “Improved interspecies optical clock comparisons through differential spectroscopy,” *Nat. Phys.* 19 (2023)

Primary author: FORTIER, Tara (NIST, Boulder)

Presenter: FORTIER, Tara (NIST, Boulder)

Session Classification: Precision Measurements and Fundamental Physics I

Track Classification: Precision Tests on Fundamental Physics