Contribution ID: 5 Type: **not specified**

Results with thin and standard p-type detectors after heavy neutron irradiation

Monday 12 November 2007 11:30 (20 minutes)

Comparison of the charge collection efficiency of $300\mu m$ and $140\mu m$ thick silicon microstrip detectors after neutron irradiation up to $1x10^{^{\circ}}16$ cm-2.

Primary author: Dr CASSE, Gianluigi (Department of Physics)

Presenter: Dr CASSE, Gianluigi (Department of Physics)

Session Classification: Full detector systems

Track Classification: Full Detector Systems