

MRI Guided Proton Therapy: Pencil beam scanning in an MRI fringe field

Brad Oborn^{1,2}, S. Dowdell³, P. E. Metcalfe^{2,4}, S. Crozier⁵, S. Guatelli²,
A. B. Rosenfeld², R. Mohan⁶, P. J. Keall^{7,4}

¹Illawarra Cancer Care Centre (ICCC), Wollongong Hospital, NSW 2500 Australia

²Centre for Medical Radiation Physics (CMRP), University of Wollongong, NSW 2522 Australia

³Shoalhaven Cancer Care Centre, Nowra, NSW 2541 Australia

⁴Ingham Institute for Applied Medical Research, Liverpool, NSW 2170 Australia

⁵School of Information Technology and Electrical Engineering, University of Queensland, St Lucia,
QLD 4072 Australia

⁶Department of Radiation Oncology, MD Anderson, Houston, Texas 77030

⁷Sydney Medical School, University of Sydney, Sydney, NSW 2006 Australia

ICTR-PHE: Geneva, Feb 15-19 2016



UNIVERSITY OF
WOLLONGONG



MRI-guided Proton Therapy (MRPT):

- A **proposed** modality that uses proton beams to treat cancer patients while being imaged inside an MRI scanner


Why consider MRPT?

- Improve the efficacy of proton beam therapy:
 - Superior 4D patient anatomy at **Planning** and **Treatment**:
 - Reduced margins, improved tumor delineation and dose calculations
- Cost effective: $\$(MRPT) \sim 1.1 \times \$(Proton\ Therapy)$

Obvious concerns with MRPT?

- **Protons will deflect in the magnetic field!**
- **MRI \iff Proton beamline interactions!**

Real-time MRI guided: promises for proton beam therapy





Adapted from Bjerre *et al.* Phys Med Biol 2013

MRPT: A prototype based on the Sydney MRI-linac



- Agilent 1 T Split-bore MRI
- Rotating couch = no proton gantry required



MRPT: background

Previous works have looked at proton beam deflection in phantoms and air:


- 90 and 90-200 MeV proton beam deflections[1,2]
- Planning study, 35 cm uniform fringe field[3,4]
- Full fringe field transport, 1 T Agilent MRI design[5]

Conclusions:



- 1. Moderate but predicatable deflection inside phatom:** 28 mm for 200 MeV at 1.5 T
- 2. Strong but predicatable deflection from fringe field:** MRI and orientation specific
- 3. Pencil Beam Scanning a must:** fringe field disperses multiple energy broad beams

1. Raaymakers *et al.*, Phys. Med. Biol. 53(20), 5615-5622 (2008).
2. Wolf and Bortfeld, Phys. Med. Biol. 57(17), N329-N337 (2012).
3. Moteabbed *et al.*, Med. Phys. 41(11), 111713 (11pp.) (2014).
4. Hartman *et al.*, Phys. Med. Biol. 60, 5955-5969 (2015).
5. Oborn *et al.*, Medical Physics, 42, 2113 (2015).


MRPT: Inline Orientation fringe field deflection





300 MeV diverging pencil beam phase space source
(emulates beam scanning)



MRPT: Perpendicular Orientation fringe field deflection



300 MeV diverging pencil beam phase space source
(emulates beam scanning)



Aim + Methods

Aims:


- To model pencil beam scanning in MRPT and assess:
 - MRI image distortion
 - Scanning magnet operation
 - MRI and Scanning Magnet cross-talk
 - Magnetic forces

Methods:

- COMSOL Finite element model: Agilent 1 T MRI + IBA Nozzle Scanning magnets in clinical position
 - Produce field maps for Monte Carlo model
 - Calculate image distortion
 - Calculate forces on scanning magnets
- Geant4 Monte Carlo model:
 - Pencil beam deflection process through combined system
 - Pencil beam dose to water phantom in magnetic field



Methods: The core finite element model

- MRI coils, pencil beam scanning magnet coils, and steel yokes



Methods: Pencil beam scanning = changing current density in scanning magnets


- Model a scanning calibration pattern:
 - 230 MeV proton beam, 81 points over a 30x30 cm field
- Model MRI off and on
- Each scan point simulation: 8 cores, 20 hours, 60 GB RAM



Results: MRI model

Excellent match to Agilent specs:

- 6.8 parts per million (ppm) magnetic field non-uniformity in the 30 cm diameter spherical volume (DSV)
- contour match at all values




Far Fringe


Near Fringe

Results: MRI + Scanning magnets model (passive)

- Inline Orientation: 6.8 ppm \mapsto 56.5 ppm
- Perpendicular Orientation: 6.8 ppm \mapsto 120.5 ppm
- Both of these are $<$ 300 ppm preshim threshold (correctable)




Inline



Perp


Results: Inline Orientation Scanning Pattern

230 MeV pencil beam scanning, scanning magnets in clinical position: (gif)





Results: Perpendicular Orientation Scanning Pattern

230 MeV pencil beam scanning, scanning magnets in clinical position: (gif)



Results: Example full beam transport to a water phantom

- To get a Bragg-peak at 190 mm depth (Perpendicular Orientation)....



- Its clear that a full trajectory model is required for planning purposes

Results: Summary



- MR image distortion (extreme case):
 - Inline Orientation: 56.5 ± 0.75 ppm
 - Perpendicular Orientation: 120.5 ± 0.4 ppm
 - Shimable, expect no image distortion
- Pencil Beam Scanning Process:
 - Scanning magnets work as intended...
 - Minimal cross-talk expected (1.5 ppm change in extreme case)
 - Forces < 1000 N
 - Fringe field corrections required for downstream deflections

Discussion and Conclusions...

- The performance of a proton pencil beam scanning assembly has been modelled in the presence of and MRI fringe field:
 - Some MRI imaging correction required, **will need active/passive shimming**
 - Scanning process appears normal, **downstream fringe correction (software/planning) required**
- Next steps including modeling detector systems

Discussion and Conclusions...

- The performance of a proton pencil beam scanning assembly has been modelled in the presence of and MRI fringe field:
 - Some MRI imaging correction required, **will need active/passive shimming**
 - Scanning process appears normal, **downstream fringe correction (software/planning) required**
- Next steps including modeling detector systems



I hope that was a real *spin!*
Are you now *distorted*
towards MRI-guided
Proton Therapy?

