

PET in Norway an example

Julie Haglund

Medical Physicist
Sweden/Norway

haglundj@aol.com

Institute For Medical Physics (IFMP) Workshop
"PET Principles, PET use in hospital, and ongoing developments"
Ohrid, Macedonia; 6-8 September 2015


Østfold Hospital Kalnes



Norway
Population: 5 165 800

Area: 323,802 km²






PET i Norge

- 3 PET/CT Oslo Universitetssykehus
- 1 PET/CT Haukeland Universitetssykehus⁺¹
- 1 PET/CT UNN
- PET/MR og PET/CT St.Olav høsten 2013

→ PET/CT Sykehuset Østfold Kalnes 2015

Østfold Hospital Kalnes

New facility
Replaces 5 regional hospitals
Psychic and somatic health care
First departments: May 2015
Complete clinical operation: November 2015



Construction start: 2011



First clinical activities:
May 2015

Insert
PET/CT
here!





What must we consider when installing a PET/CT?

- Which one to purchase?
- What will we do with it?
- What can we do with it?
- Where will it be?
- Who will use it?



Today, we have one SPECT/CT
and 2 hotlabs.
We perform nuclear
medicine studies
and diagnostic CT.



The Mo/Tc99m generator is delivered
every week.



Isotopes we use:
Tc99m
Ra-223 (therapy)
I-131 (therapy)
I-123

What do we need for PET/CT?



Positron emitting isotopes!

Our options: Positron emitting isotopes are produced in generators or cyclotrons.

*Maybe, in the future,
but also expensive.*



(Rubidium)
generator



**Out of budget!
Not spontaneously
installed!**

Cyclotron


Solution: Oslo has cyclotron

Limitations:

- What isotopes Oslo produces
- How often Oslo produces
- What can be transported (~1hr)

Conclusion:

FDG
delivery once daily



Which PET/CT to purchase?
an investment now and for the future (~10-12 years)

What can we use with a PET/CT?  FDG
Rubidium

What will we do with a PET/CT?


- Clinical patient studies/examinations
- Research projects/Collaborations

Which PET/CT to purchase?

Requires extensive evaluation
by a multi-disciplinary team!



The purchasing process includes research by a multi-disciplinary team including physicians, physicists, building/construction representatives, medical technical personnel, radiographers/technicians who will work with the PET/CT.



The group gathers information about available PET/CT systems and evaluates according to a list of desired and required characteristics of the machine that will be purchased.

Specification of Requirements

- Big Excel document
- Given to each manufacturer
- Our demands and desires for the performance of the PET/CT
- Includes:
 - physics of the PET, physics of the CT, hardware, software, IT-solutions, reconstruction algorithms, protocols, energy use, heating/cooling requirements, area required for the scanner and additional equipment, education of personnel, equipment updates,...
- Each manufacturer states how or if thier machine meets our list of requirements



Specification of Requirements

E98		fæ			
A	B	C	D	E	
109	5.33	Effektivitet (geometrisk) i prosent	B	JA	0,93
110	5.34	Detektorbredde i aksialretning (mm)	B	JA	19,2 mm.
111	5.35	Oppgraderingsmuligheter	B		Siemens har som policy at våre nyvinninger også skal kunne bli tilgjengelig på eksisterende utstyr. Det betyr at man gis mulighet til å oppgradere systemet med ny funksjonalitet om det blir aktuelt.
112	6.0	Scan			
113	6.1	Oppgi minste mulige kollimering ved aksialt og helikalt opptak	B	JA	0.6 mm Merk! Oppløsning på 0.33 mm (x,y og z) oppnås ved overlappende utregninger av snitt.
114	6.2	Oppgi største mulige kollimering ved aksialt og helikalt opptak	B	JA	20 mm
115	6.3	Oppgi maksimal mAs på lengste scan.	B	JA	666 mAs
116	6.4	Oppgi detektordekning ved helikalt og aksiale skann	B	JA	128 mm
117	6.5	Bordhøyde (mm) min/max	B	JA	Table speed range 0.1-200 mm/sec ; gjelder mCT Flow (opsjon)
118	6.6	Oppgi pitch muligheter	B	JA	0,3 - 1,5 for Biograph mCT 64,
119	6.7	Rotasjonshastighet ved 360 grader	B	JA	64 : 0.33, 0.5, 1.0 s
120	6.8	Maksimal skanntid (sek)	B	JA	Max scantid på 100s ved spiral scann i praksis ingen begrensning.
121	6.9	Maksimal skannlengde (cm), både helikalt og aksialt	B	JA	2200 mm for CT only Maximum co-scan: 203cm with pallet extensions (3-ring) 198cm with pallet extensions (4-ring)
122	6.10	Maks antall bilder mulig i en serie	B	JA	Ingen praktisk begrensning. Man kan for eksempel kvele hele kroppen, med 0.33 mm i romlig oppløsning. Etter dette kan dere rekonstruere så mange bilder i alle planer i ønsket snittykkelse, alt etter hva slags kollimering dere har valgt. Og i tillegg kommer det an på hvor stort mellomrom dere ønsker mellom snittene. Fra én rotasjon med ett rør/én detektor kan systemet beregne 384 snitt (1mm).
123	6.11	Rekonstruksjonstid (bilder per sek)	B	JA	Opptil 40 bilder/sek for Biograph mCT 64
124	6.12	Minste tidsoppløsning (ms)	B	JA	Ned til 83 ms (Biograph mCT 64)

Requirement

Importance scale

Fulfilled or not

Details: how the requirement is fulfilled

Side 1


Which PET/CT to purchase: 3 choices



GE Discovery 610



Siemens Biograph mCT



Philips Vereos PET/CT

Which PET/CT to purchase: Evaluating the PET Physics



Characteristics we wanted:

The scanner shall be able to perform all types of PET examinations and come with protocols and equipment for performing all scan types.


Evaluation:

- The machine from every manufacturer gives the possibility to perform all relevant examinations.
- Siemens & GE: requested software can be purchased as an option
- Philips: all requested software is included; protocols from our Ingenuity CT can be imported to the PET/CT

Characteristics we wanted: EKG and respiratory gating

Evaluation:

- Each manufacturer has a unique solution for EKG and respiratory gating, and the solution is included in the purchasing of the PET/CT
- Siemens: a unique system for gating and triggering that makes possible the gating of CT and PET uptake; EKG triggering is possible; gating is prospective or retrospective
- GE: uses Varian RPM and IVY7800 trigger for cardiac gating
- Philips: is compatible with Varian RPM



Evaluating the Physics

Evaluating the PET Physics means looking at the crystals, the detectors

Properties of the Ideal Scintillation Crystal for PET


Crystal property	Purpose
High density	High γ -ray detection efficiency
High atomic number	High γ -ray detection efficiency
Short decay time	Good coincidence timing
High light output	Allows large number of crystal elements per photodetector
Good energy resolution	Clear identification of full energy events
Emission wavelength near 400 nm	Good match to photomultiplier tube response
Transparent at emission wavelength	Allows light to travel unimpeded to photomultiplier tube
Index of refraction near 1.5	Good transmission of light from crystal to photomultiplier tube
Radiation hard	Stable crystal performance
Nonhygroscopic	Simplifies packaging
Rugged	Allows fabrication of smaller crystal elements
Economic growth process	Reasonable cost

Melcher, Charles L.; Scintillation Crystals for PET; J NuclMed2000;41:1051-1055

Characteristics: Scintillation crystals



Scintillator characteristics

Manufacturer	Scintillator crystal	Light decay time (ns)	Attenuation coefficient (cm ⁻¹)	Light output (photons/MeV)	Density (g/cm ³)
GE	BGO	300	0,95	7000	7,13
Siemens, Philips	LSO, LYSO	40	0,86	26000	7,4



- Crystal density influences detection sensitivity
- Light decay time influences count rate; determines scanner deadtime and random coincidences rate; determines use of ToF or not
- Light output determines energy and spatial resolution
- Stopping power of LSO is slightly lower than BGO

Evaluating the Physics



Manufacturer	Detector technology	Crystal size (mm)	# crystals / detector
GE	Analogue	4,7 x 6,3 x 30	12 crystals/PM-tube
Siemens	Analogue	4 x 4 x 20	42 crystals/PM-tube
Philips	Digital	4 x 4 x 19	1:1 coupling (23,040 digital detectors)

- Larger crystals \Rightarrow wider FWHM; worse spatial resolution
- Crystal elements coupled to PM-tubes \Rightarrow worse spatial resolution, even for a greater total number of detector elements
- Smaller crystal cross-sectional area \Rightarrow smaller pixel size; greater spatial resolution



Spatial Resolution (FWHM)

a measure of the ability of the device to faithfully reproduce the image of an object; the minimum distance between two points in an image; influenced by detector size & block configuration

Siemens

PET NEMA 2007 Spatial Resolution – Axial (Typical)	Standard Processing (256x256)	HI-REZ Processing (400x400)
FWHM @ 1 cm (mm)	5.5	4.5
FWHM @ 10 cm (mm)	6.0	5.9
PET NEMA 2007 Spatial Resolution – Transverse (Typical)	Standard Processing (256x256)	HI-REZ Processing (400x400)
FWHM @ 1 cm (mm)	5.9	4.4
FWHM @ 10 cm (mm)	6.0	4.9

Philips

Transverse spatial resolution @ 1 cm	4.1 mm FWHM
Transverse spatial resolution @ 10 cm	4.5 mm FWHM
Axial spatial resolution @ 1 cm	4.1 mm FWHM
Axial spatial resolution @ 10 cm	4.3 mm FWHM

GE

Spatial resolution FWHM ¹	NEMA performance standards ^{††}	VUE Point HD ¹
Axial @ 1 cm	5.6 mm	5.0 mm
Axial @ 10 cm	6.3 mm	5.0 mm
Transaxial @ 1 cm	5.0 mm	4.0 mm
Transaxial @ 10 cm	5.6 mm	4.5 mm




System Sensitivity

- the number of counts per unit time detected by the PET/CT for each unit of activity present in a source
- depends on the geometric efficiency, detection efficiency, PHA window settings, and the dead time of the system
- detection efficiency of a detector depends on the scintillation decay time, density, atomic number, and thickness of the detector material

Manufacturer	System Sensitivity (cps/kBq)
GE	10
Siemens	9,5
Philips	22

System sensitivity increase as the square of the detector efficiency, which depends on scintillation decay time and stopping power of the detector → LSO/LYSO preferred to BGO





Axial detector width

More detector material across the FoV captures more information per bed position

Manufacturer	Axial FoV (cm)	Number of detector rings
GE	15,7	4
Siemens	21,6	4 (TrueV)
Philips	16,4	5

We desire as wide a detector as possible in the axial direction → fewer bed stops to cover the same patient volume





Bed overlap


PET/CT signal suffers from degradation from the center to the edge of the axial FoV

- Overlapping sequential bed positions compensates
- Can lead to axially varying noise sensitivity if there is not enough overlap

Manufacturer	Bed Overlap (%)	Options, Adjustments
GE	23% (standard 11 slices)	Overlap can be adjusted with respect to 1-23 slices
Siemens	43%	FlowMotion eliminates overlap; maintains uniform noise sensitivity across the entire scan
Philips	32%	Optimized for the scanner

Time of Flight

Smaller time resolution, coincidence window, and energy resolution → better time and spatial resolution



Manufacturer	Timing Resolution (ps)	Coincidence Window (ns)
GE	550	9,5
Siemens	540	4,1
Philips	345	4,0

GE offers a reconstruction method Q.Clear to process the data instead of measuring ToF

Noise Equivalent Count Rate

Image noise is characterized by NECR

$$\text{Noise Equivalent Count Rate (NECR)} = \frac{\text{TRUES rate}^2}{\text{TRUES rate} + \text{RANDOMS rate} + \text{SCATTER rate}}$$


- NECR is proportional to SNR in the final reconstructed images
- NECR serves as a good parameter to compare the performances of different PET scanners
- Image noise can be minimized by maximizing NECR

Noise Equivalent Count Rate

SIEMENS

COUNTS VS AKTIVITETSKURVE – VEREOS PET

NEC – clinical performance




r sec
90
60

INJECTED DOSE RANGE
185 – 740 MBq
5 – 20 mCi
1 hour uptake

Biograph HI-REZ PICO

KING HAMAD
UNIVERSITY HOSPITAL
مستشفى الملك حمد الجامعي

3D NECR Curves.



Diagnostic Acquisition
NU2-2007 3D SF and Cnt Rates Scan 1
nuclideName: 18F
data taken at 13/3/2012 16:19:00
report made at 29-Mar-2012 10:04:51

peak NECR is 79.6 KHz at 12.6 kBq/cc
SF at peak NECR is 39.17 %

3D Peak NECR

	Counts [kcp/s] ^[1]	Activity Conc [kBq/cc] ^[2]
SPEC	68 kcp/s	15 kBq/cc
Measured	79.6 kcp/s	12.6 kBq/cc
Pass/Failed	Pass	N.A.

[1] - SPEC : Lower Limits
[2] - Measured concentration at peak NECR depends on activity meter calibration and frame sampling time. More than 20% variation in the reported concentration is expected. This value is not a performance specification.



- GE Confidential -

11 /
GE Bartholomé V. /
3/30/2012

Noise Equivalent Count Rate

- NECR describes the effective number of counts measured as a function of activity in the FoV
- $NECR \propto SNR$
- Image noise can be minimized by maximizing NECR

Scatter Fraction

Lower SF \Rightarrow better scanner performance
Lower SF \Rightarrow better quality of images

$$SF = \frac{C_s}{C_p}$$

Scatter, Prompt count rates

Manufacturer	Peak NECR	Transaxial FoV/bore diameter (cm)	Scatter fraction (%)
GE	76kcps@15kBq/mL	70	38%
Siemens	175kcps@ \leq 28kBq/mL	78	38%
Philips	650kcps@50kBq/mL	67,6	31%



Which PET/CT to purchase?

A	B	C	D	E
5.33	Effektivitet (geometrisk) i prosent	B	JA	0,93
5.34	Detektorbredde i aksialretning (mm)	B	JA	19,2 mm.
5.35	Oppgraderingsmuligheter	B		Siemens har som policy at våre nyvinninger også skal kunne bli tilgjengelig på eksisterende utstyr. Det betyr at man gis mulighet til å oppgradere systemet med ny funksjonalitet om det blir aktuelt.
6.0	Scan			
6.1	Oppgi minste mulige kollimering ved aksialt og helikalt opptak	B	JA	0.6 mm Merk! Oppløsning på 0.33 mm (x,y og z) oppnås ved overlappende utvalg av snitt.
6.2	Oppgi største mulige kollimering ved aksialt og helikalt opptak	B	JA	20 mm
6.3	Oppgi maksimal mÅs på lengste scan.	B	JA	666 mÅs
6.4	Oppgi detektordekning ved helikale og aksiale skann	B	JA	128 mm
6.5	Bordbevegelse (mm/s) min/max	B	JA	Table speed range 0.1-200 mm/sec ; gjelder mCT Flow (opsjon)
6.6	Oppgi pitch muligheter	B	JA	0.3 - 1.5 for Biograph mCT 64,
6.7	Rotasjonshastighet ved 360 grader	B	JA	mCT 64 : 0.33, 0.5, 1.0 s
6.8	Maksimal skanntid (sek)	B	JA	Max scantid på 100s ved spiral scann
6.9	Maksimal skannlengde (cm), både helikalt og aksialt	B	JA	I praksis ingen begrensning. 2200 mm for CT only. Maximum co-scan: 203cm with pallet extensions (3-ring) 198cm with pallet extensions (4-ring)
6.10	Maks antall bilder mulig i en serie	B	JA	Ingen praktisk begrensning. Man kan for eksempel kjøre hele kroppen, med en 1000 mm i romlig oppløsning. Etter dette kan dere rekonstruere så mange bilder dere ønsker i ønsket snittykkelse. Etter hva slags kollimering dere har valgt. Og tillegg kommer det an på hvor stort mellomrom dere ønsker mellom snittene. Fra én rotasjon med ett rør/én detektor kan systemet beregne 384 snitt (1mm).
6.11	Rekonstruksjonstid (bilder per sek)	B	JA	Opptil 40 bilder/sek for Biograph mCT 64
6.12	Minste tidsoppløsning (ms)	B	JA	Ned til 83 ms (Biograph mCT 64)



Side 14

Our team members evaluated their parts of the specification of requirements and presented for the group, and the purchasing package from every manufacturer was evaluated on an equivalent basis

Purchasing a PET/CT



**Insert
PET/CT
here!**

Philips Vereos PET/CT received the highest evaluation and was purchased.

first digital PET/CT in Norway, third Vereos installation in Europe



Purchasing a PET/CT




GE Discovery 610
11 450 000 NOK
(128slice): +400000NOK
Medrad Intego
939 000NOK



Siemens Biograph mCT 64
12 500 000NOK
Flow: 16 750 000NOK
True V: +1621250NOK
Intego: 836 600NOK


Philips Vereos PET/CT
(64slice): 17 825 000NOK
(128slice):+1200000NOK
Intego: 849 800NOK





Install Vereos...

...but first install
Gemini 64!



Gemini 64, September 2015
PET patients, October 2015
Vereos, April 2016(?)

What must we consider when installing a PET/CT?

Where will it be located?

We have a room (actually 2 rooms) shielded for SPECT/CT


Tc99m energy: 140keV

halflife: 6,02 hours


SPECT shielding is insufficient for PET/CT!

F-18 energy: 511keV

halflife: 109,8 min




A disadvantage of deciding to install PET after building for SPECT:



shielding for SPECT \neq shielding for PET
→ rooms get smaller

shielding for SPECT \neq shielding for PET


3mm lead doors into the patient injection rooms



2mm lead SPECT/CT

4mm lead doors into the PET/CT camera room

shielding for SPECT \neq shielding for PET



Østfold Hospital Kalnes has a
delay tank in the basement

Emptied every 24 hours

What must we consider when installing a PET/CT?

Who will use the PET/CT?

Our staff is well trained and routined in SPECT/CT:

myocardial scintigraphy

CBF

heart stress/rest

diagnostic CT with contrast

I-131 thyroid therapy

Ra-223 Xofigo therapy

I-123

Tc99m



All staff are new to PET/CT!

Tc99m energy: 140keV
halflife: 6,02 hours




F-18 energy: 511keV
halflife: 109,8 min

Norway is relatively new to PET/CT

- 4 PET/CT Oslo Universitetssykehus
- 1 PET/CT Haukeland Universitetssykehus
 - 1 PET/CT UNN
- PET/CT og PET/MR St. Olsv høsten 2013
- PET/CT Sykehuset Østfold Kalnes 2015

Survey March/April 2014:
How did other hospitals in Norway
educate their staff?



Survey: How did other hospitals educate their staff?



- EANM courses in Wien
- Manufacturer education upon installation
- Study visits to other hospitals
- Reading/studying articles
- Department seminars from colleagues who had been to other hospitals, physicists




Much emphasis on EANM guidelines and study visits



- We also sent out staff to other hospitals in Norway and Denmark
- Asked for protocols and patient letters used by other hospitals to be the guide for our own documents
- EANM guidelines
- Department seminars by the physicist
- Physicist went to many conferences, courses



The SPECT/CT staff knows:



What is different with PET/CT:



Learn examination protocol from EANM:

Community Medical Center
Department of Radiology

Patient History Questionnaire


Patient Name: _____ Date: _____
DOB: _____ Weight: _____

Reason for Procedure:

Please check any of the following symptoms related to today's test:

- Chest Pain Headaches Nausea Hearing Loss
 Abdominal Pain Blackouts Blurred Vision Ringing in Ears
 Pelvic Pain Dizziness Memory Loss Neck Pain
 Back Pain Unexpected Weight Loss
 Shoulder Pain (Right/ Left) Leg Pain (Right/ Left)
 Arm Pain (Right/ Left) Numbness (Right side/ Left side)
 Weakness (Right side/ Left side)

How and when did these symptoms occur (e.g., injury, just started, etc.) _____



LEFT FRONT BACK RIGHT
Please identify the location of any pain/numbness/lump

Medical History:

- Do you have or have you had any of the following?
 Cancer Kidney/Renal Disease Multiple Myeloma Tumor, lump or mass
 Bleeding Tendency Diabetes Congenital Heart Defect Stroke
 Hypertension Seizures Sickle Cell Anemia
 Asthma, bronchitis or emphysema Other illness/disease: _____
 - Have you had any tests (MRI, CT, X-Ray, etc.) performed for the symptoms you are currently experiencing? Yes No
 If yes, please list the date, type and who performed the test: _____
 - Have you had any surgeries or therapies (e.g., radiation therapy, chemotherapy, etc.) Yes No
 If yes, please list the date and type of surgery or therapy: _____
 - Are you currently taking any medications? Glucophage Metformin Glyburide
 If yes, please list all medications you are currently taking on the attached form.
 - Do you have any allergies (e.g. medications, latex, food, etc.)? Yes No
 If yes, please list all allergies: _____
 - When is your next appointment with the physician that ordered this procedure? _____
 - Females only: Are you or is it possible you may be pregnant? Yes No Don't know
 Breast Feeding? Yes No First day of last period (LMP) _____
- I hereby certify that the above information is true and to the best of my knowledge.

Patient or Legal Representative Signature _____ Print Name and Authority (if legal representative) _____ Date _____

Technologist Notes: _____

Eur J Nucl Med Mol Imaging (2010) 37:181–200
DOI 10.1007/s00259-009-1297-4

GUIDELINES

FDG PET and PET/CT: EANM procedure guidelines for tumour PET imaging: version 1.0

Ronald Boellaard · Mike J. O'Doherty · Wolfgang A. Weber · Felix M. Mottaghy · Markus N. Lonsdale · Sigrid G. Stroobants · Wim J. G. Oyen · Joerg Kotzerke · Otto S. Hoekstra · Jan Pruim · Paul K. Marsden · Klaus Tatsch · Corneline J. Hoekstra · Eric P. Visser · Bertjan Arends · Fred J. Verzijlbergen · Josee M. Zijlstra · Emile F. I. Comans · Adriaan A. Lammertsma · Anne M. Paans · Antoon T. Willemsen · Thomas Beyer · Andreas Bockisch · Cornelia Schaefer-Prokop · Dominique Delbeke · Richard P. Baum · Arturo Chiti · Bernd J. Krause

Published online: 14 November 2009

© The Author(s) 2009. This article is published with open access at Springerlink.com

Abstract The aim of this guideline is to provide a minimum standard for the acquisition and interpretation of PET and

PET/CT scans with [18F]-fluorodeoxyglucose (FDG). This guideline will therefore address general information about [18F]-fluorodeoxyglucose (FDG) positron emission tomography-computed tomography (PET/CT) and is provided to help the physician and physicist to assist to carrying out, interpret, and document quantitative FDG PET/CT examinations, but will concentrate on the optimisation of diagnostic quality and quantitative information.

This guideline is a joint project of the EANM Oncology Committee and the EANM Physics Committee. In addition, this guideline is based on the following three documents:

(1) DGN (Deutsche Gesellschaft für Nuklearmedizin) Leitlinie: "FDG-PET/CT in der Onkologie" by Krause BJ, Beyer T, Bockisch A, Delbeke D, Kotzerke J, Minkov V, Reiser M, Willich N, Arbeitsausschuss Positronenemissionstomographie der Deutschen Gesellschaft für Nuklearmedizin, 2007.

(2) SNM Guidelines: "Procedure Guidelines for tumour imaging with 18F-FDG PET/CT 1.0." by Delbeke D, Coleman RE, Guiberteau MF, Brown ML, Royal HD, Siegel BA, Townsend DW, Berland LL, Parker JA, Hubner K, Stabin MG, Zubal G, Kachelries M, Cronin V, Holbrook S, 2006.

(3) "Applications of F18-FDG-PET in Oncology and Standardisation for Multi-Centre Studies" by Boellaard R, Oyen WJG, Hoekstra CJ, Hoekstra OS, Visser EP, Willemsen AT, Arends AJ, Verzijlbergen JF, Paans AM, Comans EFL, Lugtenburg E, Stoker J, Schaefer-Prokop C, Zijlstra JM, Pruim J. HOVON Imaging workgroup and the Netherlands Society of Nuclear Medicine, 2007

R. Boellaard (✉) · O. S. Hoekstra · E. F. I. Comans ·

Keywords Guideline · FDG · PET · PET/CT · Tumour · Oncology · Quantification · QC · QA

Introduction

The aim of this guideline is to provide a minimum standard for the acquisition and interpretation of PET

W. A. Weber

Create a patient questionnaire based on examples from other hospitals

EANM EARL Accreditation



Is your
PET/CT centre
accredited
by the EANM
Research Ltd.
(EARL)?

Østfold Hospital Kalnes:

- regional hospital, not a university hospital
- new to PET/CT

EARL accreditation:

- demonstrates that we know what we are doing
- possibly facilitates research collaborations
- third EARL hospital in Norway

Means more work for the physicist!


- 4 quality control tests (image quality, recovery coefficients) to submit per year

Takk!



Julie Haglund
Medical Physicist
haglundj@aol.com


Characteristics: scan range and timing window for reconstruction



Manufacturer	Shortest frame-time	Maximum number of frames	Scan range
GE	1 second	250	120cm
Siemens	3 seconds	100	195cm with TrueV and pallet extensions
Philips	6 seconds randoms subtraction included	99	190cm



<http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3826954/>



The central image compares TOF and non-TOF PET imaging. On the left, a circular detector ring is shown with a grid of detector elements. A central inset shows a PET scan of a human torso. On the right, a similar setup is shown with a stopwatch displaying '495 ps', representing the time resolution of the TOF system. The text below explains that TOF technology uses the actual time difference between the detection of coincident events to more accurately identify the origin of the annihilation, leading to a quantifiable improvement in image quality.

Time
Time

TOF

No TOF

TOF

495 ps

Time-of-Flight technology uses the actual time difference between the detection of coincident events to more accurately identify the origin of the annihilation. Better identification leads to a quantifiable improvement in image quality.