



Canadian Association
of Physicists

Association canadienne
des physiciens et physiciennes

Contribution ID: 2134

Type: **Oral (Graduate Student) / Orale (Étudiant(e) du 2e ou 3e cycle)**

Optical fiber nanoprobe for the detection of chemicals (G)

Wednesday, 13 June 2018 14:45 (15 minutes)

Surface Enhanced Raman Spectroscopy (SERS) is a powerful tool for detecting chemicals, including rhodamine 6G (R6G), at the molecular level. The authors have developed a plasmonic structure by depositing gold nanorods (GNRs, aspect ratio: 6.4) on a tapered optical fiber via optical tweezing. A 1064 nm laser was coupled to the un-tapered end of the tapered fiber for tweezing. The tweezing process was repeated with several other laser wavelengths such as He-Ne 632 nm. In order to show the efficiency of the plasmonic structure the authors will present the SERS spectra of R6G. Authors will also describe techniques (Transmission, and reflection modes) which can be used in conjunction with the nano-probe in order to collect scattered light from the sample.

The research was financially supported by Natural Sciences and Engineering Research Council of Canada (NSERC) and Canada Foundations for Innovations. Authors also acknowledge the members of the Photonics Research Group, Dr. Reznik's research group, and Dr. Aicheng Chen's research group.

Primary authors: Dr DAS, Gautam (Lakehead University); Mr TREVISANUTTO, Joshua (Lakehead University)

Presenter: Mr TREVISANUTTO, Joshua (Lakehead University)

Session Classification: W3-5 Quantum Sensors and Quantum Technology (DAMOPEC) | Détecteurs quantiques et technologie quantique (DPAMPC)

Track Classification: Division of Atomic, Molecular and Optical Physics, Canada / Division de la physique atomique, moléculaire et photonique, Canada (DAMOPEC-DPAMPC)