



Contribution ID: 1129

Type: **Parallel Session Talk**

Measurement of the top quark mass and width with CDF detector

Friday 23 July 2010 11:25 (22 minutes)

We report a measurement of the top quark mass obtained from proton-antiproton collisions at a center-of-mass energy of 1.96 TeV at the Fermilab Tevatron using the CDF II detector using different decay modes and analysis methods. We present measurements of the top quark width, and mass difference between the top and anti-top quark in the lepton+jets channel as well.

Author: CDF COLLABORATION, CDF (Fermi National Accelerator Laboratory (FNAL))

Presenter: LEE, Hyunsu (University of Chicago)

Session Classification: 02 - The Standard Model and Electroweak Symmetry Breaking

Track Classification: 02 - The Standard Model and Electroweak Symmetry Breaking